

# Active II Water Birth Pool





Judges at the Building Better Healthcare awards said that the original Active Birth Pool was an incomparable combination of intelligent design, superior materials and state-of-the-art manufacturing.

The latest advancements in design, material and engineering take the new model of this award winning water birth pool to new levels of safety, usability and end-user satisfaction.



## Specifications

Length:	1965mm
Width:	1650mm
Height:	770mm
Depth:	640mm
Weight (empty):	100kg
Capacity:	560 litres
Material:	Ficore® composite
Guarantee:	Lifetime



# Active II Water Birth Pool

ACTIVE  
BIRTHPOOLS

Multi-colour Lighting and Blue-tooth Sound create a conducive atmosphere for physiological labour and natural birth

The organic shape of this luxurious pool is appealing both physically and psychologically

Seamless one-piece fabrication in Ficore composite denies bacteria a place to grow

The distinctive keyhole shape channel protects mothers from hitting their heads on the spout

Recessed hand grips do not obstruct movement around the rim and provide support in upright postures

Superior material highly resistant to bacteria and chemicals

The highly polished surface is ultra-smooth, less slippery and warm to the touch

The wide wrap-around rim maximises comfort and usability for mothers and midwives

Portable hoist compatible

There is maximum space under the rim for midwives to sit closely in comfort

The bull-nose rim provides a firm hand hold for mothers in all positions

Simple and safe for mothers to get in and out of the pool

The Safety Seat facilitates monitoring and emergency evacuations

Optimum space, and depth and give mothers complete freedom of movement

The Labour Support Seat provides support in upright positions and a place for mothers to rest and welcome their baby

