

ELITE

4 CHANNELS ELECTROTHERAPY | 98 PROGRAMS

FOUR CHANNELS PRODUCT FOR EVERY TREATMENT NEEDS.

32 SPORT | 29 BEAUTY | 26 FITNESS | 9 TENS | 2 REHAB

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

4 independent (8 electrodes)

Frequency

0,3-150Hz *

Pulse amplitude

40-450 μ s *

Power

0-100mA per channel

Power supply

Rechargeable batteries

Size

mm 100x160x35


Weight

gr 454

* According to the type of current

EQUIPMENT

- 1 Bag
- 1 Stimulator ELITE
- 4 Cables for electrode connection
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE
0476



TENS



LAST 10



MULTI USER



MY TRAINER



SYNCR0 STIM



RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

ELITE

4 CHANNELS ELECTROTHERAPY | 98 PROGRAMS

PROGRAMS LIST

SPORT

Demo

Motor point pen

Capillarization

Warm-up

Decontracting

Post-competition recovery/training

Maximum strength

Endurance strength

Explosive strenght

Reactivity

FITNESS

Firming

Bio-Pulse Firming

Sculpting

Toning

Mass building

BEAUTY

Drainage

Bio-Pulse Drainage

Lipolysis

Bio-Pulse Lipolysis

Toning massage

Bio-Pulse relaxing massage

Vascularization

Face Capillarization

Face lifting effect

TENS AND PAIN

Antalgic tens

Endorphinic tens

Contractures

Hematomas

Sciatica

Lumbago

Cervical pain

Shoulder pain

Epicondylitis

REHAB

Lower limbs atrophy recovery

Upper limbs atrophy recovery

