

# PHYSIOLASER 1000

## LASER THERAPY

PHYSIOLASER 1000 IS A PRODUCT OF THE PROFESSIONAL LINE. POWERFUL AND PORTABLE, THE PRODUCT IS INDICATED TO ACCELERATE THE RECOVERY OF INFLAMED TISSUES, REGENERATE CELLS AND IMPROVE MICROCIRCULATION.

### TECHNICAL SPECIFICATIONS

#### Display

2,6", backlit

#### Laser source

Class 4 GaAlAs diode

#### Maximum power

1W  $\pm$ 20%

#### Compatible with a wavelength of

808nm

#### Issue

Pulsed and continuous

#### Target area dimensions

11 mm optic

#### Frequency

Up to 10,000Hz

#### Power supply


Mains and rechargeable batteries

#### Size and weight

100x160x35 mm - 433 g

### STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 1000 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



CE  
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

**GLOBUS**  
ITALIAN EXCELLENCE

# PHYSIOLASER 1000

## LASER THERAPY

### PROGRAM LIST

#### HEALTH - REHAB

Arthrosis

---

Ankle sprain

---

Knee pain

---

Ankle sprain edema

---

Meniscopathy

---

Knee osteoarthritis

---

Achilles tendinitis

---

Finger pain

---

Myofascial shoulder pain

---

Articular shoulder pain

---

Epicondylitis

---

Frozen shoulder

---

Carpal tunnel

---

Chronic pain Back pain

---

Acute back pain

---

Back pain with radiculopathies

---

Neck pain

---

Neck pain with radiculopathies

---

Myofascial neck pain

---

Neck osteoarthritis

---

Temporo-mandibular pain

---

Tendinopathies

---

Injuries

#### AESTHETIC - FITNESS

##### FACE AESTHETIC

Wrinkles

---

Acne

---

Dehydrated skin

---

Skin rejuvenation

---

Skin - face lift

---

Expression wrinkles

---

Neck wrinkles

---

Rosacea

---

Telangiectasias

##### BODY AESTHETIC

Breast stretch marks

---

Edematous cellulite on thighs, hips and gluteus

---

Fibrous cellulite on thighs, hips and gluteus

---

Fat and non-localized cellulite

---

Stretch marks

---

Epilation

---

Scars

---

Loc. lymphatic stasis