



# THE BOOST 2 CORE

International Technical Specifications

- Speed 0-19 kmh
- Incline 0-15%
- Weight-bearing 100-20%
- Color White
- User height 1.20m - 2.05m
- Footprint 2.08m L x 1.06m W
- Power Requirements 120V/20A or 220V/20A
- Touchscreen display 58cm; shows % weight reduction, walking speed, race pace per km, distance, heart rate
- Data Storage Via the Cloud (no USB needed)
- Gait Analysis Via Digitsole, clinical gait views
- CE & FDA Approval CE approved & FDA cleared
- Shorts 2XS - 5XL, compatible with AlterG®
- User controls Simplified, ergonomic knobs
- Base treadmill Woodway 4Front Treadmill
- Height Adjustment Automatic & precise based on user height
- Workouts Pre-programmed or design your own

