

# PHYSIOLASER 12.0 PRO

## LASER THERAPY


THE "PRO" VERSION OF PHYSIOLASER 12.0 IS DESIGNED FOR THE PROFESSIONAL WHO USES THE DEVICE FREQUENTLY AND PREFERS TO KEEP IT ALWAYS NEARBY ON A CART. THE 12W OF MAXIMUM POWER ALLOW RAPID AND EFFECTIVE TREATMENTS EVEN ON THE LARGEST BODY AREAS, USING THE 22 MM OPTIC INCLUDED. THE WIDE LIBRARY OF PROGRAMS ALSO GUARANTEES THE PROFESSIONAL TO TREAT ANY MUSCULOSKELETAL DISORDER AND SOME PROBLEMS RELATING TO WOUNDS AND ULCERS OF THE SKIN.

### TECHNICAL SPECIFICATIONS

<b>Display</b>
2.6", backlit
<b>Laser source</b>
Class 4 GaAlAs diode
<b>Maximum power</b>
12W ±20%
<b>Compatible with a wavelength of</b>
808nm 980nm and 1064nm (optional)
<b>Issue</b>
Pulsed and continuous
<b>Target area dimensions</b>
11 mm optic 22 mm optic
<b>Frequency</b>
Up to 10,000Hz
<b>Power supply</b>
Mains and rechargeable batteries
<b>Size and weight</b>
220x170x60 mm - 790 g


### STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 12.0 PRO device
- 1 808nm laser handpiece with 11mm and 22mm optics
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



### PURCHASE OPTIONS

- Code G5995 HP 980nm handpiece
- Code G5996 1064nm handpiece, maximum power 3W



# PHYSIOLASER 12.0 PRO

## LASER THERAPY

### MEDICAL PROGRAM LIST

#### 808nm WAVELENGTH PROGRAMS

##### MUSCULOSKELETAL

###### GENERAL

Arthrosis

Rheumatoid arthritis

Muscle injuries

Muscle contracture

Decontracting

Fibromyalgia

Radiculopathies

Tendinopathies

Neuralgia

Postherpetic neuralgia

Sudeck disease

Bursitis

Edema

###### ARM - HAND

Arthrosis

Rheumatoid arthritis

Raynaud syndrome

Finger pain

Epicondylitis

Carpal tunnel

Wrist fractures

Wrist pain

Tendinitis

Edema

Radius fracture

###### KNEE

Arthrosis

Rheumatoid arthritis

Acute pain

Chronic pain

Meniscus pain

Bursitis

Patello-femoral

Edema

Ligament injury

Gonalgia

Sprain

Chondromalacia

Chondropathy

Patello-femoral tendinitis

###### FOOT - ANKLE

Arthrosis

Rheumatoid arthritis

Metatarsalgia

Ankle pain

Plantar fasciitis

Achilles tendon

Sprain

Edema

Heel pain

###### BACK

Vertebral arthrosis

Acute pain

Lumbago

Myofascial pain

Radiculopathies

Rheumatoid arthritis

###### SHOULDER

Myofascial pain

Arthrosis

Adhesive capsulitis

Bursitis

Subacromial impingement

Tendinitis

Luxations

Tendon injuries (cuff)

Rheumatoid arthritis

Edema

###### HEAD - NECK

Arthrosis

Radiculopathies

Myofascial pain

Chronic cervical

Temporo-mandibular pain

Rheumatoid arthritis

Trigeminal neuralgia

Cervical herniation

##### DERMATOLOGY

Ulcers

Cicatriziation ulcers

Wounds

#### 980nm WAVELENGTH PROGRAMS

Antalgic

Cervical Pain

Lumbar pain

Bursitis

Baker's cyst

Cruralgia

De Quervain syndrome

Sprains

Contusions

Acute pain

Chronic pain

Edemas

Hematomas

Extended hematomas

Epicondylitis

Epitrochleitis

Discal hernia

Fibromyalgia

Inflammations

Ligament injury

Muscle injuries

Meniscus pain

Synovitis

Muscle strains

Tendinopathy

Carpal Tunnel