

# CEZANNE JR 3

Physiotherapy table with three sections.



Medical  
Device



Dimensions  
197 x 64 x 50 h92 cm



Adjustable height  
da cm 50 cm min a 92 cm max



Safe working load  
160 Kg



Maximum capacity  
200 Kg



Unimpeded passage  
for patient lifting



The photo depicts the product complete with optional accessories

## FEATURES

---

Electric motor with a thrust force of 8000 N

---

Push-button control panel

---

---

Class 1IM fire-retardant padding and upholstery

---