

Nearly 25 years ago, NuStep created the first product of its kind – the recumbent cross trainer. Now a mainstay in health, wellness, and fitness settings, the recumbent cross trainer has transformed hundreds of thousands of lives worldwide.

NuStep recumbent cross trainers offer an ideal workout solution for users who are fit, deconditioned, living with chronic conditions or recovering from a health event.

Today, NuStep continues its mission to transform lives with the launch of T6. Built on our proven T5 cross trainer platform, the T6 delivers the same cardio and strength workout users rely on, but with advanced features and programs to enhance their workout experience.

Why NuStep?

Quality:

We manufacture NuStep cross trainers at our headquarters in Ann Arbor, Michigan USA to ensure they meet strict quality standards.

Reputation:

Our superior products and longevity in the marketplace have earned us a reputation for excellence among industry professionals in healthcare, senior care and fitness.

Satisfaction:

Our number one goal is customer satisfaction. We focus on the quality and durability of our cross trainers to ensure they withstand years of demanding use.



+1 734.769.3939 | nustep.com NuStep®, StrideLock®, and WellGrip® are registered trademarks of NuStep, LLC, Ann Arbor, Michigan, USA • 2020
The Bluetooth® word mark logo is a registered trademark owned by Bluetooth SIG, Inc. and any use of such marks by NuStep, LLC is under license.
Other trademarks and trade names are those of their respective owners.



T6: STEP BEYOND EVERY EXPECTATION.




Features

Every detail was considered.

- > 15 levels of resistance, easily changed during your workout using the lower screen
- > Ergonomic hand grips rotate 40 degrees and sturdy arm handles, with easy-to-use clamshell release levers, adjust up to 18" (46 cm)
- > 12-degree seatback recline
- > Durable, oversized foot pedals have non-slip inserts and adjustable foot straps
- > Integrated media shelf
- > Compatible with T5 series accessories and adaptive equipment
- > Easy start up – system wakes up by touching the screen or stepping
- > Water bottle holder included
- > Optional grab ring and hip and torso belts for additional stability

A smart, sensible interface.

- > The upper screen provides an excellent interface for user engagement with an enhanced Pace Partner program, Quick Start and Charting Apps
- > The lower screen displays all workout metrics including time, speed (steps per minute), total distance, watts, calories, METs, resistance level and heart rate*
- > Admin mode allows for customization of the default settings
- > Heart rate monitoring through **Bluetooth®** wireless technology*
- > Audible feedback with volume control
- > 3-minute console shutdown with extended pause option for interval training
- > Seat position readout displayed on screen
- > USB port for data export to PDF



StrideLock® feature locks handles and pedals



12-degree seatback recline offers maximum comfort



Adjustable foot straps to fit every user

Specs

Machine footprint: 73"L x 30"W x 52"H (185 cm x 76 cm x 117 cm)
Weight: PRO: 285 lbs. (129 kg) / MAX: 298 lbs. (135 kg)
Standard seat width: 18" (46 cm)
Wide seat width: 22" (56 cm)
Fits users from 4'6" to 6'7" (137 cm to 200 cm) in height
Power: Medical grade 24V DC power supply included

Build your T6

Base Model Part No.	Description
06001	T6 _{PRO} with standard 18" (46 cm) width seat
06003	T6 _{MAX} with 22" (56 cm) width seat

Optional Add-ons*

Part Number	Description
52341	Grab ring support
52421	Hip belt
59095	Torso belt

Adaptive Equipment

Part Number	Description
52000	Unilateral Leg Stabilizer Kit Includes 1 leg stabilizer, 1 strap (26"/66 cm) and 2 storage clips
52065	Bilateral Leg Stabilizer Kit Includes 2 leg stabilizers, 2 leg straps (26"/66 cm) and 2 storage clips
50680	WellGrip® – left hand only
50681	WellGrip – right hand only
50685	WellGrip Set – left and right hand

Accessories

Part Number	Description
50780	Accessory Storage Bag
41815	Floor Mat
05300	NuStep Transporter
51877	Polar®H7 Heart Rate Sensor & Strap

*Must be ordered at time of T6 purchase. Not available as a retrofit option.

*When a Bluetooth heart rate monitor is connected.