

Xcite2

Task-Specific Activities to Optimize Functional Goals

Xcite2

Advanced Technology for Coordinated iFES™ Evoked Movement Patterns

Xcite2 was designed to allow clinicians to easily and efficiently apply multi-channel FES with task-specific activities commonly completed in therapy.

Xcite2 has 40+ pre-programmed activities including:

- Activities of Daily Living (ADL) skills
- Pre-gait skills
- Mobility skills
- Fine motor skills

Individuals with neurological impairments have the opportunity for quality, task-specific repetitions of meaningful functional activities to optimize recovery potential. Xcite2 supports clinicians and patients in this important part of ongoing rehabilitation.



restorative-therapies.com



Xcite2 is Indicated for General Rehabilitation for:

- Muscle re-education
- Relaxation of muscle spasms
- Prevention or retardation of disuse atrophy
- Increasing/maintaining local blood circulation
- Maintaining or increasing range of motion

Technical Specifications:

- Height of cart: 29"–45"
- Weight of cart: 32.5lbs
- Weight of Xcite2: 6lbs
- Usable work space: 30"
- Tablet: 10.5" Touch Screen

 Restorative
THERAPIES
Your advocate for life.

Xcite2 Technology with Optimization



Features:

Xcite2 Clinical Station provides a variety of coordinated, iFES™ evoked, dynamic movement patterns available in four activity libraries:

- Hand activities
- Upper extremity activities
- Lower extremity activities
- General activities

Xcite2 Activities Include:

- Grasp & release
- Feeding
- Transfers
- Side-stepping
- Sit to stand
- And many more functional skills

- 40+ pre-programmed iFES™ functional activities
- Patient settings sync with RTILink.com
- Upgraded controller provides increased processing speed and improved touch interface
- Removable tablet allows for flexibility in treatment options
- Long-lasting battery life for 8 hours of continuous use between charges
- Height-adjustable cart allows for standing or sitting positions
- Large lockable wheels for easy maneuvering

Restorative Therapies Inc.

8098 Sandpiper Circle, Suite M, Nottingham, MD 21236

restorative-therapies.com

+1 (800) 609-9166 | Fax: +1 410 878 2466

To locate a distributor, please visit restorative-therapies.com/international.

If your country is not listed, please contact support@restorative-therapies.com.

Xcite2 is for clinical use only. Services offered vary by country.

Contact your local distributor for complete offerings in your area.

Xcite2 Sales Sheet English V0



Restorative
THERAPIES

Your advocate for life.



CE
1639